

2235D

Facial Chair / Podiary Chair



SILVERFOX
beauty

Silverfox Corporation Limited

TEL: 0086-750-2633320 2633805

E-mail: sales@silverfox.cn

ADD: No.18, 1st TongLe Road, TangXia Town,
PengJiang District, JiangMen City, GuangDong Province, China

© silverfoxbeauty1995 Silverfox Beauty

CE FDA RoHS

www.silverfoxgp.com

UNPARALLELED USER EXPERIENCE



※ Standard Configuration

- Max load bearing
- Reset

※ Options


- Heating
- Latex / Memory Foam
- PVC bed Cover
- Memory System
- Hand Controller

※ Technical Parameters

Model	Product	Voltage	Frequency	Power Rating	Warranty
2235D	Facial Chair Podiary Chair	AC 100-120V	50/60Hz	300W	2 years (Motors and Control Box: 3 years)

※ General Specifications

- Max load bearing : 175 kg/386 lbs
- Net weight : 106.5kg/ 235 lbs
- Gross weight : 141 kg/ 310 lbs
- Product size(L*W*H) : 1840x850x590/880 (mm) 72.4x33.5x23.2/34.6 (in)
- Package size : 1480x720x640 (mm)≈0.68m³ 81.1x33.9x24 (in)≈0.68m³
- Unit:mm / in



LEATHER COLOR SWATCH

STANDARD LEATHER



MICROGROOVE LEATHER



MICROFIBER LEATHER (OPTION)




NEW LEATHER (OPTION)



Maintenance

After every treatment, please only use neutral solution to clear the dust, stains on leather, metal part or other parts. Disinfecting solution like alcohol, acid or alkaline solution are forbidden. This is to prevent the lifetime of the components from being shortened.

SPONGE DISSECTION DIAGRAM



Due to the screen display effect will have color difference, the image is for reference only, the actual color shall prevail.